FREE ACTION

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GUIDANCE AND MANUFACTURER'S DECLARATION ELECTROMAGNETIC EMISSIONS

The Free Action is intended for use in the electromagnetic environment specified below. The customer or the user of the Free Action should assure that it is used in such an environment.

EMISSIONS	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
RF emissions CISPR 11	Group 1	The Free Action uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Free Action is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Conform Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Conform	

GUIDANCE AND MANUFACTURER'S DECLARATION ELECTROMAGNETIC IMMUNITY

The Free Action is intended for use in the electromagnetic environment specified below. The customer or the user of the Free Action should assure that it is used in such an environment.

IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines and dc lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
	±1 kV for input/output lines	Not applicable for input/output lines	
Surge IEC 61000-4-5	±1 kV line(s) and neutral ± 2 kV line(s) to earth	±1 kV line(s) and neutral N/A(57.2 EN61000-4-5)	Mains power quality should be that of a typical commercial or hospital environment. If a dips or an interruption of mains power occurs, the current of the Free Action may be dropped off from normal level, it may be necessary to use uninterruptible power supply or a battery.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % Ur (>95 % dip in Ur) for 0,5 cycle	<5 % Ur (>95 % dip in Ur) for 0,5 cycle	
	40 % Ur (60 % dip in Ur) for 5 cycles	40 % Ur (60 % dip in Ur) for 5 cycles	
	70 % Ur (30 % dip in Ur) for 25 cycles	70 % Ur (30 % dip in Ur) for 25 cycles	
	<5 % Ur (>95 % dip in Ur) for 5s	<5 % Ur (>95 % dip in Ur) for 5s	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE - Ur is the a.c. mains voltage prior to application of the test level

GUIDANCE AND MANUFACTURER'S DECLARATION ELECTROMAGNETIC IMMUNITY

The Free Action is intended for use in the electromagnetic environment specified below. The customer or the user of the Free Action should assure that it is used in such an environment.

IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT – GUIDANCE
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3V	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Free Action, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p style="text-align: center;">Recommended separation distance $d = 1.17 \sqrt{P}$ $d = 1.17 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.33 \sqrt{P}$ 80 MHz to 2.5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation Distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	Interference may occur in the vicinity of equipment marked with the following symbol:

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations.
Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Free Action is used exceeds the applicable RF compliance level above, the Free Action should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Free Action.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE FREE ACTION

The Free Action is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Free Action can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Free Action as recommended below, according to the maximum output power of the communications equipment.

RATED MAXIMUM OUTPUT POWER OF TRANSMITTER (W)	SEPARATION DISTANCE ACCORDING TO FREQUENCY OF TRANSMITTER (M)		
	150 KHZ TO 80 MHz $d = 1.17 \sqrt{P}$	80 MHz TO 800 MHz $d = 1.17 \sqrt{P}$	800 MHz TO 2,5 GHZ $d = 2.33 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 - At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 - These guidelines may not apply in all situations.
Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

FREE ACTION

CARTE DE GARANTIE / WARRANTY CARD / GARANTÍA / GARANTIE-KARTE 24 MOIS / 24 MONTHS / 24 MONATE / 24 MESES / 24 MAANDEN

Nom / Name / Name / Nombre /
 Nome / Naam / Nome / Jméno /
 Név / nazwa / nume / име /
 namn / имя / isim / اسم

Adresse / Address / Anschrift / Dirección /
 Indirizzo / Adres / Endereço / Adresa /
 Cím / adres / adresa / адрес /
 address / адрес / adres / عنوان

Code Postal / Postcode / Postleitzahl / Código Postal /
 Codice postale / Postcode / CEP / PSČ /
 Irányítószám / Kod pocztowy / Cod poştal / Пощенски код /
 Postnummer / почтовый индекс / posta kodu / الرمز البريدي

Ville / Place / Stadt / Ciudad /
 Città / Stad / Cidade / Město /
 Város / miasto / oraş / град /
 Stad / город / kasaba / بلدة

Date d'achat / Date of purchase / Einkaufsdatum / Fecha de compra /
 Data di acquisto / Verkoopsdatum / Data da compra / Datum nákupu /
 A vásárlás dátuma / Data zakupu / Data de cumpărare / Дата на закупуване /
 Inköpsdatum / Дата покупки / Satin alma tarihi / تاريخ الشراء

CACHET DU REVENDEUR / STAMP OF RETAILER / STEMPEL VOM HÄNDLER / SELLO DEL VENDEDOR /
 TIMBRO DEL RIVENDITORE / STEMPEL VAN DE VERKOPER / CARIMBO DO REVENDEDOR / RAZÍTKO PRODEJCE /
 STAMP A KISKERESKEDELMI / PIECZĘĆ SKLEP / STAMPILA DISTRIBUTOR / ПЕЧАТ НА МАГАЗИНИ /
 STÄMPEL ÄTERFÖRSÄLJARE / ПЕЧАТЬ ДИСТРИБЬЮТОРЫ / SATINÇININ STAMP / الطابع من متاجر التجزئة

En cas de panne retournez l'appareil accompagné de cette carte de garantie à votre revendeur.
 In case of damage, please return appliance and warranty card to your retailer.
 Im Garantiefall, senden Sie bitte das Gerät mit der Garantiekarte an unsere Kundendienststelle.
 En caso de defectos, por favor entregue al vendedor el equipo junto con esta garantía.
 En caso di danni, restituire l'apparecchio con questa carta di garanzia al rivenditore.
 In geval van defect stuur het artikel terug naar de verkoper vergezeld van deze garantiekart.
 Em caso de danos, o retorno do aparelho com este cartão de garantia para o revendedor.
 V prípade porúbení, vráťte prístroj s tímto záručným listu na svého predajcu.
 Abban az esetben, kár, vissza a készüléket a garancialevél a kiskereskedő.
 W przypadku szkody, zwróć urządzenie z kartą gwarancyjną do sprzedawcy.
 În caz de avarie, întoarceți aparatului cu acest card de garanție pentru distribuitorul dumneavoastră.
 В случае на повреда, върнете уреда с тази гаранционна карта да си търговец на дребно.
 I händelse av skada, returnera apparaten med denna garanti kort till din återförsäljare.
 В случае повреждения, вернуть устройство с этим гарантийным талон к вашему продавцу.
 Hasar durumunda, perakendedi bu garanti karti ile cihaz dönün.
 في حالة وقوع ضرر ، وعودة هذا الجهاز مع بطاقة الضمان لتاجر التجزئة.

SPORT-ELEC® INSTITUT

B.P. 35 - 31 rue du Val Breton - 27520 BOURGTHEROULDE - FRANCE
 Tel. +33 (0)2 32 96 50 50 - Fax. +33 (0)2 32 96 50 59 - info@sport-elec.com

IP22

CE
0459

Made in P.R.C.

SPORT-ELEC® INSTITUT

B.P. 35 - 31 rue du Val Breton - 27520 BOURGTHEROULDE - FRANCE
Tel. +33 (0)2 32 96 50 50 - Fax. +33 (0)2 32 96 50 59 - info@sport-elec.com